International Instructions for Use

XLight

Languages:

EN, DE, NL,
NO

fEX Engineering GmbH

Höltengeweg 49
48155 Münster, Germany

Fourth edition

Publication date 22/07/2024

Foreword

Thank you for purchasing an **XLight** lighting system from fEX Engineering GmbH.

The XLight is an accessory for the Power-PRO XT wheeled stretcher made by *Stryker*. The system consists of a right and a left unit. It helps emergency services personnel when working at night by providing bright surround lighting, and its flashing warning light makes the scene of an incident safer.

WARNING

Improper handling or any unintended use may cause incur risk or cause damage. Please therefore read and carefully follow these Instructions for use. Retain them for reference at all times. To prevent personal injury or damage to property, please observe all safety precautions.

To contact the manufacturer:

fEX Engineering GmbH
Höltenweg 49
48155 Münster, Germany
engineering@fex.group
+49 (0)251 508 537 -13

Contents

1	How to use these Instructions for use.....	1
1.1	Feedback on these Instructions for use	1
1.2	Explanation of the warning notices used in these Instructions for use	1
1.3	Abbreviations	2
2	Intended purpose	3
2.1	Intended use.....	3
2.2	Intended medical purpose	3
2.3	Clinical benefits	3
2.4	Indications	3
2.5	Contraindications	3
2.6	Intended target patient group	3
2.7	Intended use with other products.....	3
2.8	Intended users.....	4
2.9	Intended environment for use	4
3	Important safety notice.....	5
3.1	Battery	5
3.2	What to do in the event of damage	6
3.3	Electromagnetic compatibility (EMC).....	6
3.4	Reporting obligation.....	6
4	Product description	7
4.1	Product overview	7
4.2	Symbols used.....	8
4.3	How it works.....	9
4.4	Controls	9
4.5	Use in combination with other products	10
4.6	Delivery.....	10
4.7	Modifications.....	10
4.8	Requirements for storage, handling and environmental conditions	10
4.9	Product service life	10
4.10	Warranty.....	10

5	Assembly and putting into service	11
5.1	Requirements for operators.....	11
5.2	General	11
5.3	Assembly and installation.....	11
5.4	To check for correct assembly.....	12
6	Using the product.....	13
6.1	Operator’s duty of care	13
6.2	User’s duty of care.....	13
6.3	Instructions for safe operation.....	13
6.4	Introduction.....	13
6.5	How to switch on.....	14
6.6	Mode selection	14
6.7	To switch off:.....	14
6.8	Automatic power-off.....	14
6.9	Process overview	14
6.10	Charging the battery.....	15
7	Reprocessing: Cleaning and disinfection.....	16
7.1	Cleaning the product	16
7.2	Disinfectants.....	16
7.3	Disinfecting the product:.....	17
7.4	Final treatment.....	17
8	Maintenance.....	18
8.1	Maintenance.....	18
8.2	Inspection	18
8.3	Troubleshooting	19
8.4	Damage.....	20
8.5	Maintenance and repairs	21
9	Dismantling, recycling, disposal	22
9.1	Dismantling.....	22
9.2	Recycling, disposal.....	22
9.3	Disposal through fEX Engineering GmbH	22
10	Risk management.....	23
10.1	Foreseeable misuse	23
10.2	Implemented risk control measures	23
10.3	Residual risk.....	23

11	Performance data.....	24
11.1	Product conformity	24
11.2	Manufacturer’s declaration of electromagnetic immunity.....	25
12	Technical data.....	26
12.1	Left unit and right unit.....	26
12.2	Battery	26
12.3	LED per unit	26
12.4	Connection and charging cable	26
12.5	Charger	27

1 How to use these Instructions for use

If you have any questions about the contents of these Instructions for use or how to use the product, please contact us. You can request a paper copy of these instructions for use or download the current, electronic version from our website:

<https://www.fex.group/downloads>

Should you need further assistance after carefully reading these Instructions for use and additional information, please contact your authorised distributor or contact us directly.

1.1 Feedback on these Instructions for use

Your opinion is important to us. Please let us know if there is anything you would change about these instructions for use (engineering@fex.group). We will take your feedback into account in the next version of these Instructions for use.

1.2 Explanation of the warning notices used in these Instructions for use

Signal word	Description
 WARNING	Potentially hazardous situation resulting in severe injury or death.
 CAUTION	Potential hazard resulting in minor or moderate injury.
NOTE	User error can damage the product.

Table 1: Explanation of the warning notices used in these Instructions for use

1.3 Abbreviations

Abbreviation	Definition
BfArM	Bundesinstitut für Arzneimittel und Medizinprodukte (Federal Institute for Drugs and Medical Devices)
BMS	Battery management system
EMC	Electromagnetic compatibility
ESD	Electrostatic discharge
HF communication devices	High-frequency communication devices
IFU	Instructions for use
IP	International protection
LED	Light-emitting diode
ME devices/ME systems	Medical electrical devices/medical electrical systems
MPSV	Medizinprodukte-Sicherheitsplanverordnung (Medical Devices Safety Plan Ordinance)
UDI	Unique Device Identification

Table 2: Abbreviations

2 Intended purpose

2.1 Intended use

The XLight is designed as an extension of the *Power-PRO XT* stretcher made by *Stryker* to light the surrounding area and make the point of use safe with a warning light.

Hand-held use away from the wheeled stretcher is prohibited.

2.2 Intended medical purpose

The XLight is an active, non-invasive Class I medical device.

2.3 Clinical benefits

The XLight is classed as an accessory and has no inherent clinical benefit. It is designed to help the *Power-PRO XT* wheeled stretcher made by *Stryker* to fulfil its intended purpose.

2.4 Indications

No indications have been specified for the XLight.

2.5 Contraindications

There are no known contraindications for the XLight.

2.6 Intended target patient group

The patient population is actually using the *Power-PRO XT* wheeled stretcher made by *Stryker*. No direct contact with the patient is anticipated.

2.7 Intended use with other products

The XLight is intended exclusively as an accessory for the *Power-PRO XT* medical device made by *Stryker*.

2.8 Intended users

The XLight is intended for use in a professional healthcare environment. The users, and the qualifications they require, are specified accordingly.

User groups	Emergency services personnel (primary users)	Medical technicians, technical managers, cleaning staff (secondary users)
Training, certification	Medical training, e.g. as an emergency paramedic, emergency medical technician, rescuer worker, emergency services support staff	Medical technicians (training/academic qualifications) or equivalent, or accreditation as a cleaner in a professional healthcare environment

Table 3: Intended users

2.9 Intended environment for use

The XLight is intended for use by primary users in the following environments:

- Use in proximity to patients and therefore in a professional healthcare environment
- Use in a mobile environment, e.g. in emergency services vehicles
- Use outdoors
- Temperatures of -20°C to 50°C
- IP65 (dust proof and waterproof)

3 Important safety notice

NOTE

- Please read and observe the instructions in these Instructions to use this system safely.
- Do not open the product.

⚠ CAUTION

- Only use it with the connection cable provided by the manufacturer. Otherwise, EMC performance may be impaired and interference with other ME equipment and ME systems may occur.
- Keep away from extreme heat sources, such as open flames. Using the device within a temperature range of -20°C to +50°C.
- The XLight is not approved for use in the air service.

⚠ WARNING

- The XLight must be fully locked into place. Otherwise there is a risk that the XLight could detach from the *Power-PRO XT* in the event of a traffic accident and become a projectile that could injure the occupants.
- The XLight is not intended for use in an oxygen-rich environment.

3.1 Battery

The XLight contains a rechargeable LiFePO₄ battery (12 V/3 Ah) made by *Just Battery Technology UG* with a built-in battery management system (BMS). You can find technical details about this battery in Section 12.2 *Battery*. The batteries have been tested by the manufacturer in accordance with UN38.3.

Only the charger from fEX Engineering may be used to charge the battery. You can find technical details about the charger in Section 12.5 *Charger*.

⚠ CAUTION

All batteries carry the risk of fire, explosion or leakage of hydrofluoric acid. Possible reasons:

- Extreme heat sources, such as open flames
- Temperature range beyond -20°C to +50°C
- Severe impacts or crushing
- Short circuit

⚠ CAUTION

Only the manufacturer may open this product. This also applies to replacing the battery.

3.2 What to do in the event of damage

The XLight must be withdrawn from use if it is damaged (Section 8.4 *Damage*). Repairing damaged components is prohibited – these must be replaced (Section 8.5 *Maintenance*).

Visible damage must be immediately reported to the technical manager.

3.3 Electromagnetic compatibility (EMC)

Due to its electrical components, the XLight emits electromagnetic radiation which can interfere with other ME systems.

 **CAUTION**

Only use it with the connection cable provided by the manufacturer.

 **WARNING**

The XLight must be mounted in the location provided for it (Section 5.3 *Assembly and mounting*). Other ME systems must not be mounted or carried within 30 cm.

Failure to comply with these warnings may impair EMC performance and endanger safety. For more information see Section 13 *Manufacturer's declaration - Electromagnetic immunity*.

3.4 Reporting obligation

All serious incidents which occur in connection with the product must be reported to the manufacturer and the competent authority in the EU Member State where you are established.

In Germany, the competent authority is BfArM. Up-to-date contact details can be found on the BfArM website: <https://www.bfarm.de>

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4 Product description

4.1 Product overview



	<p>XLight</p> <ul style="list-style-type: none"> • Patient right unit • Both units are symmetrical <ol style="list-style-type: none"> 1. 3 LED 2. 3 reflective rings 3. Control button 4. Locking unit
	<p>Connection cable</p> <ul style="list-style-type: none"> • Angled plug on both sides • Both plugs are identical • 0.4 m in length
	<p>Locking unit</p> <ol style="list-style-type: none"> 1. Positioning aid 2. Twist knob 3. Batch number 4. Lock position (<i>open/closed</i>)
	<p>Charger and -charging cable</p> <ul style="list-style-type: none"> • Charger for simultaneous battery charging <ol style="list-style-type: none"> 1. Mains adaptor 2. Control LED 3. Y-charging cable 4. Power cable

Table 4: Product overview

4.2 Symbols used


Figure 1: Right unit cover with safety information

Symbol	Meaning
	Manufacturer's logo
	Manufacturer, address and contact details
Charging input	Charging with 14.6 V and 2 A
	DC
	Warning
	Medical device
	Conformity marking
1.4 kg	Weight per lighting unit
	Read Instructions for use
IP 65	Protection class
	Proper disposal of old electrical and electronic equipment
	Battery supply with 12.8 V direct current
	Permissible temperature range -20°C to +50°C
	UDI – Unique Device Identification: (00) UDI-DI (Device Identifier; (10) Batch no.; (11) Date of manufacture

Table 5: Symbols used and their meanings

4.3 How it works

The XLight is an accessory for the *Power-PRO XT* wheeled stretcher made by *Stryker*, and consists of one right unit and one left unit. It helps emergency services personnel when working at night by providing bright surround lighting, and its built-in warning light warning light makes the scene of an incidentsafer.

Both units are fit with positive locking onto the wheeled stretcher plate and locked underneath with quick-release fasteners.

Each unit has three high-power LEDs pointing in different directions:

- Direction of travel – white (in patient’s line of sight)
- Surroundings – white (almost 360° when both units are in use)
- Warning light – yellow/white (flashing)

Controlling the three modes is intuitive because the buttons on both units are easy to reach and synchronised with each other.

- Mode 1: standard surround lighting (direction of travel + surroundings)
- Mode 2: maximum surround lighting (direction of travel + surroundings)
- Mode 3: Warning light (constant in direction of travel plus yellow/white flashing light)

Both sides have independent controls and their own battery. They each weigh 1.4 kg. The 3000 mAh LiFePO₄ battery can provide continuous lighting for up to 24 hours.

The system meets the standards of protection class IP65, is approved under the current Medical Devices Ordinance (MDR) and has been successfully dynamically tested according to standard DIN 1789 (10g/5 axes).

4.4 Controls

Control	Function	Operation	Location
	On and off button Selects the mode	To switch on: Press the switch To select mode: Press the switch To switch off: Hold the switch down for 3 seconds	See Figure 2: XLight product overview – patient right, no. 4 (page 9)
	Locks the units onto the wheeled stretchers	To close: Turn clockwise to “ON” To open: Turn anticlockwise to “OFF”/🔓	See Figure 3: Component view of locking unit, no. 2 (page 9)

Table 6: Controls

4.5 Use in combination with other products

The XLight is designed exclusively for use with the *Power-PRO XT* wheeled stretcher made by *Stryker*. The product can be installed and locked onto this wheeled stretcher with positive locking. The conditions for adequately safe use are only guaranteed when used in this way.

4.6 Delivery

The product is delivered in non-sterile, non-disinfected condition and not fully charged.

An electronic copy of the Instructions for use can be found via the following link:

<https://www.fex.group/downloads/>

A sticker on the system shows a QR code and the above link, both of which will take you to the electronic version of the Instructions for use. This sticker must be removed before use.

4.7 Modifications

Unauthorised modifications to this product are prohibited.

4.8 Requirements for storage, handling and environmental conditions

Category	Information
Permitted humidity limits	5% to 98%
Permitted temperature limits (use)	-20°C to +50°C
Permitted temperature limits (storage)	-20°C to +40°C
Permitted temperature limits (charging)	0°C to +40°C
Protect the product from impacts, crushing or rough handling	

Table 7: Storage during transport and before and after use

4.9 Product service life

Under normal usage conditions and with appropriate maintenance, the service life is **seven years**.

The date for the product’s service life is based on the delivery date from fEX Engineering.

4.10 Warranty

The XLight is subject to a statutory warranty of **two years** under EU and German law.

The battery is excluded from this warranty.

5 Assembly and putting into service

5.1 Requirements for operators

The operator and his personnel can assemble the product by following the instructions below. It is the operator's responsibility to have qualified personnel carry this out.

5.2 General

The battery is delivered uncharged or partially charged and must be fully charged before first use. For instructions on how to charge the battery, see Section 7.8 *Charging the battery*.

The XLight is delivered in non-disinfected condition, and must therefore be disinfected before first use. For disinfection instructions, see Section 8 *Reprocessing: cleaning and disinfection*.

5.3 Assembly and installation

The XLight is designed exclusively for use with the *Power-PRO XT* wheeled stretcher made by *Stryker*. The positive locking mechanism on the product means it can only be used safely with this wheeled stretcher.

The lighting system consists of two units, patient left and patient right. These are connected via a connection cable. Both units are attached at the foot-end under the patient lying surface in the bracket weld (Figures 2 and 3) where the stretcher lock is located.

The LED system should point in the direction of the patient's line of sight. Positive locking prevents incorrect positioning/fitting.


Figure 2: Fitting position

Figure 3: Bracket weld

Figure 4: After fitting

The explanation of how to fit the unit applies to patient side right. Fitting is the same for the left side.


Figure 5: Locking unit


Figure 6: Fitting the locking unit


Figure 3: Inserting the connection cable

Follow the steps below:

1. Turn the twist knob (2) anticlockwise to "OFF / " to unlock the fastening on the locking unit (1).
2. Remove the locking unit (1) from the clamping plate (4) of the lighting unit.
3. Position the lighting unit in the bracket weld. The raised areas of the housing base (5) fit into the bracket weld with positive locking.
4. Reinsert the locking unit (1) into the clamping plate (4). The positioning aids (3) must line up with one other.
5. Turn the twist knob (2) clockwise to the "ON" position to lock.
6. Repeat steps 1 to 4 as required for the other side.
7. Insert the connection cable (6) into the sockets (7) of both units.
The connection cable (6) must be parallel to the edge of the housing (8) to avoid damage to the socket (7).

5.4 To check for correct assembly

Use the check list below to make sure the system is assembled or installed correctly. If any of these points do not apply, the product must be reassembled or reinstalled.

- The connection cable must be completely inserted to ensure a good connection.
- If the lock cannot be removed without turning the twist knob, the system is installed correctly.

Once assembled or installed properly, the XLight is ready for immediate use

6 Using the product

6.1 Operator's duty of care

To ensure safe operation, the operator must meet the following obligations:

- Personnel must be trained to use the XLight installed
- The Instructions for use must be read and followed
- Products must be assembled or installed properly
- Only suitable, qualified, trained specialist personnel may use the product
- All serious incidents (damage, injuries, etc.) must be reported to the manufacturer and the competent authorities

6.2 User's duty of care

To ensure safe operation, the user must meet the following obligations:

- The Instructions for use must be read and followed
- The product must be used as intended
- Damage must be reported to the technical manager

6.3 Instructions for safe operation

- Ensure product is fully locked
- Only use the product for its intended purpose
- Only use the product within the permitted temperature limits
- Only use in the intended position
- Follow EMC instructions
- Only use the manufacturer's connection cable

6.4 Introduction

The XLight has one button for each lighting unit. These buttons control the entire lighting system. Once the connection cable is connected, both batteries and both buttons are connected and synchronised. This means that the button on either side also controls the other lamp units. This applies to all of the user's actions.

The position of this button is shown Section 4.1 *Product overview*.

When the product is switched off, it goes into standby mode and discharges slowly by itself.

6.5 How to switch on

Push either of the two buttons to switch on the XLight.

Mode 1: Surround lighting with brightness dimmed

6.6 Mode selection

Press the button repeatedly to change mode.

Mode 2: Surround lighting (full brightness)

Mode 3: Warning light (surround lighting dimmed, coloured LED flashes as a warning light)

Mode 1: Surround lighting (dimmed)

The modes are programmed on a loop so the user can toggle continuously through them.

6.7 To switch off:

Press and hold one of the buttons for 3 seconds to switch off both units.

6.8 Automatic power-off

Both units automatically power off 20 minutes after the last button press.

6.9 Process overview


Figure 8: Process diagram

6.10 Charging the battery

The battery must be charged **outside** of the vehicle. We recommend charging the battery at least every 2 weeks at the same time as cleaning and inspecting the XLight. If possible, combine with existing weekly routines (e.g. when cleaning the vehicle).

The charging unit consists of a mains adapter, socket and a separate charger cable. The charging cable is designed as a Y-cable to allow both XLight units to be charged simultaneously. You can find technical details about the charger in Section 12.5 *Charger*.

To charge the battery:

1. Remove both lighting units from the wheeled stretcher.
2. Remove the connection cable from both units.
3. Plug in the charging device using the charging cable.
4. Plug the charging cable into the socket on the top of each unit.
5. Control LED on the charger:
 - a) Red: Batteries charging
 - b) Flashing red/green: Batteries are balancing
 - c) Green: Charging is complete
6. Reinstall the system, lock into place properly and reconnect the units with the connection cable.

NOTE

Make sure the charger is connected to the power supply. The built-in batteries mean the control LED lights up even when the charger is **not** connected to mains power.

WARNING

Only use the charger provided by fEX Engineering. Using other chargers may cause hazards or damage.

NOTE

Only charge the battery when the ambient temperature is between 0°C and +40°C.

7 Reprocessing: Cleaning and disinfection

The XLight is a reusable medical device. Cleaning and disinfection must follow the process in the Instructions for use, but must also be carried out after every potential contamination, however at least once a month.

We recommend cleaning and inspecting the XLight at the same time. If possible, combine with existing weekly routines (e.g. when cleaning the vehicle).

NOTE

The XLight must **not** be cleaned with high pressure. Remove the XLight when cleaning the wheeled stretcher with high pressure.

CAUTION

Clean the product with the connection cable plugged in or the contacts inside the charger socket will be exposed. Danger of short-circuit.

7.1 Cleaning the product

1. Remove the system from the wheeled stretcher.
2. Wipe the product's surfaces with a solution of mild soap and water
3. Dry all product surfaces thoroughly.

7.2 Disinfectants

Agents with the following ingredients can be used for disinfection:

- Quaternary cleaners containing less than 3% glycol ether
- Chlorinated bleaching agents (5.25 % – less than 1 part bleaching agent to 100 parts of water)
- Phenolic cleaners (active ingredient: o-phenylphenol)
- ≤ 21% isopropanol

Disinfectants made from quaternary ammonium compounds must be used in accordance with the recommended instructions for use to avoid damaging the plastic housing of the device. Quaternary cleaners and disinfectants must be used with extreme care on plastic surfaces.

WARNING

Disinfectants containing one of the following basic ingredients are not approved as they can damage the system:

- Hydrogen peroxide
- (Per-/Chlor-) acetic acid

7.3 Disinfecting the product:

1. Remove any visible contaminants from the product (Section 7.1 *Cleaning the product*).
2. Apply the recommended disinfectant solution. Do not soak.
3. Allow the product to air dry. Follow the instructions and exposure time as directed by the manufacturer of the chemical products.

Using an unsuitable agent, an excessive amount of an agent, an exposure times that is too long or failing to neutralise an agent sufficiently may damage the equipment. All information provided by the manufacturer of the agent used is decisive and must be observed at all times.

7.4 Final treatment

Residues of cleaning agents must be removed. After disinfection, always wipe the product with a clean, water-soaked cloth and allow it to air dry.

 **WARNING**

Residues of cleaning agents are corrosive and can injure users and damage the product.

8 Maintenance

The operator and his personnel can carry out the following tasks if the instructions are followed. It is the operator's responsibility to employ suitable, qualified and trained personnel carry this out.

8.1 Maintenance

The XLight is maintenance-free.

Carry out an inspection once a week (Section 8.2 *Inspection*).

8.2 Inspection

We recommend cleaning and inspecting the XLight at the same time. If possible, combine with existing routines (e.g. when cleaning the vehicle).

Functional testing and a visual inspection for damage (Section 9.3 *Damage*) must be carried out weekly.

Remove the XLight from the wheeled stretcher and check for potential damage:

- a) Functionality (lighting, locking)
- b) Screws fit flush
- c) Material damage (sharp edges, brittleness or cracks in the housing)
→ *Pay particular attention to the lock*
- d) Safety precautions no longer legible
- e) Reflective rings are missing or damaged
- f) Damage to protective insulation on the connection cable
- g) Other damage to the product

 **WARNING**

The XLight must be fully locked into place. Otherwise there is a risk of it becoming a projectile in the event of a traffic accident.

Check the charger and its cables for potential damage:

- a) Damage to protective insulation on the charger or charging cable
- b) Housing open
- c) Other damage

8.3 Troubleshooting

Users or technical managers can carry out troubleshooting. Use Table 8: *Troubleshooting*. If the problem cannot be solved or a repair is needed (Section 9.5 *Maintenance and repairs*) please contact us.

Problem	Possible cause	Troubleshooting
Product does not light up	Battery empty	Charge the battery
	Product is faulty	Return unit to manufacturer for repair
	Charging cable or charger faulty	Order replacement part (fEX Engineering GmbH) and replace components
Unit is not synchronised (control and/or flash frequency)	Connection cable not correctly inserted	Remove connection cable on one side, switch off units and check plug and socket
	First check while light is switched off, then check with connection cable connected	Remove connection cable on one side, switch off units and check plug and socket
	Connection cable faulty	Order replacement part (fEX Engineering GmbH) and replace components
Lock will not slot into place	Twist knob not in position	Turn twist knob to OFF to fit on the to ON to lock
	Locking unit faulty	Order replacement part (fEX Engineering GmbH) and replace components
	Clamping plate faulty	Return unit to manufacturer for repair
Lighting time lasts less than 24 hours	Batteries not fully charged	Charge the batteries
	Battery capacity is reduced	Return unit to manufacturer for replacement
Lighting not at full brightness/warning light does not flash	Operating temperature exceeded	Allow system to cool down
	One or more LEDs faulty	Return unit to manufacturer for repair
Battery does not charge	Charging cable or charger faulty	Order replacement part (fEX Engineering GmbH) and replace components
	Battery faulty	Return unit to manufacturer for repair

Table 8: Troubleshooting

8.4 Damage

CAUTION

Damage to the system can create a safety hazard during use and must be reported immediately to the technical manager.

The XLight must be withdrawn from use in the case of the following damage:

- Device does not start
- Battery no longer charges
- LED faulty
- Material damage (sharp edges, brittleness or cracks in the plastic housing)
- Safety precautions no longer legible
- Damage to protective insulation on the connection cable
- Other damage

WARNING

If the wheeled stretcher is withdrawn from use because of an accident, the lighting units fitted it must **not** be reused.

The charger must also be withdrawn from use if the following damage occurs:

- Damage to protective insulation on the charger
- Housing open
- Other damage

This does not apply to visual defects such as superficial scratches that cannot cause injury.

8.5 Maintenance and repairs

The following components can be replaced by the user or operator with original spare parts. These can be ordered upon request. Making repairs to the product is prohibited.

- Locking unit (Item no. fEX2000VE)
- Reflective ring (Item no. fEX2000RR)
- Connection cable (Item no. fEX2000VK)
- Charger (Item no. fEX2000LG)
- Charging cable (Item no. fEX2000LK)

 **WARNING**

Only replacement parts made by fEX Engineering GmbH may be used.
Using other parts may be dangerous.
Making repairs to the product is prohibited.

If the lighting unit is faulty or damaged, repairs may only be carried out by fEX Engineering. After contacting the manufacturer, the entire unit must be returned for repair.

The Instructions for use are decisive for all work on the product. Once damaged components have been replaced, the system can be put back into use.

fEX Engineering GmbH refuses to supply any replacement parts, carry out repairs or replace batteries for products that have reached the end of their service life.

9 Dismantling, recycling, disposal

9.1 Dismantling

 **CAUTION**

The product must be cleaned and disinfected before being dismantled. Components that have been dismantled or damaged may have sharp edges.

The operator and his personnel may dismantle the product. It is the operator's responsibility to have qualified personnel carry this out. The product is dismantled in the reverse order of its assembly.

9.2 Recycling, disposal

When disposing of the product, always follow **current local** recommendations, and/or statutory provisions on environmental protection and risk relating to recycling or the disposal of devices at the end of their useful life.

This is an electrical or electronic device that must not be disposed of in residual waste.

9.3 Disposal through fEX Engineering GmbH

If a product must be withdrawn from use due to damage or other reasons, the XLight (like all fEX Engineering GmbH products) can be returned to the company for free disposal.

10 Risk management

10.1 Foreseeable misuse

CAUTION

The system must only be used for its intended purpose. The following foreseeable misuses are prohibited:

- Use without locking or without full locking
- Use on another product than the Stryker *Power-PRO XT*
- Incorrect cleaning or disinfection
- The user changing the battery
- Modification or repair measures by the user that are contrary to these Instructions for use
- Failure to carry out regular inspections
- Use as a hand-held lighting system
- Charging from an unsuitable power source
- Use of spare parts made by another manufacturer
- Cleaning with high-pressure
- Other uses that deviate from the Instructions for use

Any improper use renders any liability or warranty from the manufacturer null and void.

10.2 Implemented risk control measures

As part of risk analysis, risk control measures to minimise risk have been defined and implemented, in accordance with the Risk Management Plan.

As far as possible, risk control measures are based on compliance with the applicable norms and standards. Technical control measures and monitoring measures have also been taken, and their effectiveness successfully verified.

10.3 Residual risk

Table 10 shows other factors to consider in order to minimise residual risk.

Residual risk	Measures to minimise residual risk
Objects become projectiles inside an emergency services vehicle	<ul style="list-style-type: none"> • Always ensure that both lighting units are fully locked. • Only use replacement parts from fEX Engineering GmbH
Battery fire or explosion -	<ul style="list-style-type: none"> • Only use the charger from fEX Engineering • Use of the connection cable is compulsory • Observe storage and usage temperatures • Only allow fEX Engineering to replace the battery
If the battery is defective, hydrofluoric acid can form	<ul style="list-style-type: none"> • Only use the charger from fEX Engineering • Only use replacement parts from fEX Engineering. • Observe storage and usage temperatures • Protect the product from severe impacts

Table 9: Measures to minimise residual risk

11 Performance data

11.1 Product conformity

The product meets the following regulatory requirements:

- (EU) 2017/745 (MDR) REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 5 April 2017 on Medical devices
- ISO 14971: 2022-04 Medical devices - Application of risk management to medical devices
- IEC 62366-1: 2021-08 Medical devices – Part 1: Application of usability engineering to medical devices
- DIN EN 1789: 2020-12 Medical vehicles and their equipment - Road ambulances
- DIN EN ISO 13485:2021-12 Medical devices – Quality management systems – Requirements for regulatory purposes
- DIN EN 60601-1: 2022-11 Medical electrical devices – Part 1: General safety requirements, including essential performance characteristics; Subsection 15.3.2 Push test
- DIN EN 60601-1: 2022-11 Medical electrical devices – Part 1: General safety requirements, including essential performance characteristics; Subsection 15.3.3 Impact test
- DIN EN 60601-1: 2022-11 Medical electrical devices – Part 1: General safety requirements, including essential performance characteristics; Subsection 15.3.6 Moulding stress relief test
- DIN EN 60601-1-2: 2022-01 Medical electrical devices – Part 1- 2: General safety requirements, including essential performance characteristics – Collateral standard: Electrical interference – requirements and tests
 - Interference emitted:
radiated interference according to DIN/EN 55011; Class B; 30 MHz – 1 GHz
 - Electromagnetic immunity:
Electrostatic discharge according to EN/IEC 61000-4-2
Radiated electromagnetic immunity according to EN/IEC 61000-4-3
Magnetic fields 61000-4-8 (sinusoidal)

11.2 Manufacturer’s declaration of electromagnetic immunity

Due to its electrical components, the XLight emits electromagnetic radiation which may interfere with other ME systems. To minimise this risk, the product and connecting cable have been tested according to DIN EN 60601-1-2. The threshold values for medical devices were met.

Guidelines and manufacturer’s declaration – electromagnetic immunity			
The XLight is intended for use in one of the electromagnetic environments below. The user must ensure that it is operated on one of these environments.			
Electromagnetic immunity test	IEC 60601 Test level	Compatibility level	Electromagnetic environment – guideline
Electromagnetic immunity test against high-frequency electromagnetic fields (according to EN 61000-4-3)	3 V/m 80 MHz - 2700 MHz 80% amplitude modulation at 1 kHz	3 V/m	Portable and HF communication devices, e.g. radios and associated accessories such as antennas can affect the XLight and should be used at a distance of at least 30 cm (12”) from the XLight and its connection cables.
Near fields of wireless RF communication devices (according to EN 61000-4-3)	27 V/m 385 MHz Pulse modulation: 18 Hz	27 V/m	Interference is possible in the vicinity of devices marked with the symbols below. 
	28 V/m 450 MHz FM+ 5 Hz Variance: 1 kHz Sine wave	28 V/m	
	9 V/m 710 MHz, 745 MHz, 780 MHz Pulse modulation: 217 Hz	9 V/m	
	28 V/m 810 MHz, 870 MHz, 930 MHz Pulse modulation: 18 Hz	28 V/m	
	28 V/m 1720 MHz, 1845 MHz, 1970 MHz Pulse modulation: 217 Hz	28 V/m	
	28 V/m 2450 MHz Pulse modulation: 217 Hz	28 V/m	
	9 V/m 5240 MHz, 5500 MHz, 5785 MHz Pulse modulation: 217 Hz	9 V/m	
Note: These guidelines may not apply in all cases. Buildings, objects and people can affect the magnitude of electromagnetic propagation.			

Table 10: Electromagnetic immunity – guidelines and manufacturer’s declaration

Failure to comply with these warnings may impair EMC performance and endanger safety. This can lead to the following risks:

- Damage to the system or other ME devices/ME systems
- One or both units may switch on or off independently or unexpectedly
- Units may switch between modes independently or unexpectedly, e.g. from surround lighting to flashing warning light

12 Technical data

12.1 Left unit and right unit

Category	Specification (per unit)
Dimensions (height x width x depth)	23 mm x 19 mm x 67 mm
Weight	1.4 kg
Housing material	PC-ABS
Button material	PA
Cover/lock material	Stainless steel
3x reflective rings	3M reflective film (RA2/C)

12.2 Battery

Category	Specification
Manufacturer	Just Battery Technology UG
Description	LiFePO ₄ battery with battery management system (BMS)
Voltage/capacity	12.8 V / 3 Ah
Number of parts	2 (1 per lighting unit)

12.3 LED per unit

Category	Specification
Manufacturer	NICHIA
2x white	297 lm, 5000 K, CRI 80
1X yellow	124 lm, amber coloured

12.4 Connection and charging cable

Category	Specification
Manufacturer	Lutronic
Connection cable	To connect both units
Charging cable	Y-cable to connect the charger to both units

12.5 Charger

Category	Specification
Mains adaptor	Just Battery Technology UG
Description	Intelligent CCCV charger for 12-V-LiFePO ₄ batteries with an end-of-charge voltage of 14.6 V
Input voltage	100-240 V _{AC} , max. 2 A, 50/60 Hz
Output voltage (too lights)	14.6 V _{DC} , 2 A

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