
Xlevator – Quick Reference Guide

⚠ WARNING

This quick reference guide does not replace the Instructions for use. It is essential to read and follow the complete Instructions for use.

1. Magnets
2. Chair plate
3. Locking button
4. Handle
5. Handle bar (*Xpedition*)
6. Bar clasp
7. Strap
8. Connection bar (*Xpedition*)
9. Wheel mount

Adjusting the height and locking

NOTE

The oil brake controls the speed of upward travel. Raise the system slowly to move it with ease.

⚠ CAUTION

The Xlevator must never be used for passenger transport. Never load the Xlevator with someone's full body weight or when it is locked.

1. Push the locking button.
2. The lock is now open and the height can be set.
3. Move the chair plate to the required height. From an ergonomic point of view, it is better to use the connecting tube on the track drive as a lifting handle.
4. Continue moving the chair plate until the lock clicks into place at the next possible point.

Empty trip

⚠ CAUTION

Hold the handle firmly with both hands to avoid a crush hazard.

1. Hold the handle with both hands and press the locking button.
2. Use your body weight to pull the system down until it locks in the middle position. Do not pull on the bar clasp of the Xlevator.
3. Put your foot on one of the wheel mounts.
4. Press the locking button and move the system to the lowest position using your body weight.

NOTE

The tension/traction/driving power gas spring continuously pulls the system upwards. Counteract this force with your body weight before opening the lock.

Removal

1. Now move the chair plate to its lowest position.
2. Loosen the straps, attach the right strap to the right magnet and the left strap to the left magnet.
3. Lift the handle bar of the *Xpedition* and roll the stair chair out of the Xlevator.
4. If necessary, the Xlevator can be raised to its highest position without the stair chair.

NOTE

The Xlevator must be in the upper position before the vehicle or door is moved. If the system is in the lowest position, the vehicle height can only be minimally lowered. This may occur when loading/entering the vehicle, or when using air suspension. Otherwise the system may be damaged.

Stowing

WARNING

The *Xpedition* must be properly stowed and fully strapped in. Otherwise there is a risk of the *Xpedition* detaching from the Xlevator and becoming a projectile in the event of a traffic accident.

1. If necessary, move the Xlevator down to its lowest position.
2. Fold the *Xpedition* and move the wheels so they are in front of the wheel mounts. The handle bar of the *Xpedition* must remain extended.
3. First lift the wheels of the *Xpedition* into the wheel mounts of the Xlevator. To use as little force as possible, hold the handle bar of the *Xpedition* with both hands and lift the *Xpedition* with your foot.
4. Both wheels must be fully seated in the wheel mounts.
5. Pull in the handle bar so that it sits behind the bar clasp.
6. Strap in the *Xpedition*. The straps must cross diagonally over the handle bar.
7. Tighten the straps firmly at the buckles.
8. Then move the Xlevator to its highest position.

Inspection

Functional tests and visual inspections for potential damage must be carried out weekly.

CAUTION

Any damage must be immediately reported to the technical manager.

WARNING


No damage to the straps must be visible, and they must be fully functional. Otherwise, there is a risk of the *Xpedition* becoming a projectile in the event of a traffic accident.

Check the Xlevator for the following functionality and possible damage:

- a) Damage to the straps and functionality of their locks
- b) Gas spring supports upward movement
- c) Oil brake slows upward movement (empty trip)
- d) Gas spring and oil brake are not leaking
- e) Unhindered linear movement (system can be moved up and down)
- f) Locking mechanism works in all positions
- g) Safety notices are legible
- h) Rust formation or other material damage (sharp edges, brittleness or cracks)
- i) Tight fit and possible damage to the holding plate
- j) Other damage to the product

Emergency mode

If any of the components fail, e.g. the locking mechanism or the linear movement, the Xlevator can be used in emergency mode to allow it to be used as a static bracket. This ensures that the vehicle can be used until repairs can be carried out. To do this, an M6x60 screw must be screwed in from the outside of the guide rails on both sides.

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